

PoolMate

User Guide — quick bench reference for choosing the correct workflow and understanding the required fields.

1. Calculate Molarities

Uses only the Sample section to convert concentration plus fragment size into nM. It does not calculate pooling volumes. **Minimum:** sample concentration + size, or direct nM.

2. Quick Pool

Creates a simple manual pooling table using total fmol or fmol per sample. It does not require platform presets, racks, buffer positions, or epMotion worklists. **Use for:** fast bench calculations and manual pooling.

3. Generate Worklists

Full epMotion workflow. Uses sample molarities, the final pool target, layout, deck settings, buffer, and tips, and exports DNA.csv and Buffer.csv. **Use for:** robot execution.

Sample section vs. Create-pool parameters

Sample section is what you currently have: sample names, wells, concentrations, sizes, and direct nM values. **Create-pool parameters** are what final pool you want: target nM or fmol, final volume, and layout.

Standalone dilution calculator

Dilute one sample or a whole list to a target concentration (nM or ng/μL). Every parameter — target, final volume, min transfer, max well/tube capacity, and rounding — has a shared default that can be overridden per row. Serial steps appear automatically when one dilution can't reach the target within the well limit.

Output files

DNA.csv: DNA transfer worklist. **Buffer.csv:** buffer transfer worklist. **Worksheet:** calculation record. **Inputs:** saved run settings.

Important: The built-in platform presets are vendor starting points. Confirm all volumes and loading concentrations against your own kit / instrument SOP before running anything on the robot.